



Mountain Rescue and Medisol

Exclusive offer

The products listed below are the most commonly ordered products. If you require any other products or accessories, do not hesitate to contact us on 01223 790124 or e-mail us at info@aedexpert.co.uk

We offer free shipping on all orders. Please note all prices exclude VAT and there is a £17 customs clearance fee per order.



Medisol
Saving Lives

Mountain Rescue packages

Rescue Package 1

DefiSign Pocket Plus Semi-Automatic with bag, rescue kit, ILCOR AED sticker and silver maintenance plan.

£1999
(excl. VAT)



Rescue Package 2

DefiSign Life Semi-Automatic with cabinet, paediatric electrodes, rescue kit, ILCOR AED sticker and silver maintenance plan.

£999
(excl. VAT)



AED's



The **DefiSign Pocket Plus Defibrillator Semi-Automatic** is a compact and versatile defibrillator that offers a range of useful features. Its **small size** makes it incredibly convenient to transport, especially when using the included bag. The defibrillator is designed to cater to users of all skill levels, making it suitable for various situations. The LCD screen of the defibrillator has an automatic light adjustment feature, ensuring optimal visibility in any environment. This screen displays clear visual instructions to guide users through the resuscitation process. An **ECG-display** option is available, which provides additional insight into the heart rhythm during and after the resuscitation process. This feature makes the defibrillator particularly well-suited for use by skilled rescuers.

The **DefiSign Life Semi-Automatic** defibrillator is perfectly designed for use by everyone. This competitively priced device gives clear verbal instructions, talking the rescuer through each step of the resuscitation process. The defibrillator starts working as soon as the front cover is opened, with clear spoken instructions, supported by images that are highlighted by LED lights. The defib explains what actions need to be taken, step-by-step, making the DefiSign LIFE exceptionally **user-friendly**.

- 10 year warranty
- Multi-lingual capabilities
- High dust & water resistance



Extras

Silver Maintenance Plan



Together we make sure that your device is always ready to use. The **Silver defibrillator maintenance plan** ensures a professional, annual check of your AED. Your device will be serviced in our **ISO13485 certified Medisol Service Centre** - by trained professionals. All transport costs are included.

The continuous service includes:

- AED status overview in service portal
- Free electrodes/battery after deployment
- 10% discount on spare parts

The professional check of the Silver AED maintenance plan includes:

- ✓ Functional test
- ✓ Shock test
- ✓ Replacement of spare parts
- ✓ Error correction
- ✓ Software updates
- ✓ Cleaning of device
- ✓ Visual inspection
- ✓ Detailed maintenance report
- ✓ Data storage for 5 years
- ✓ Maintenance advice

£119
(excl. VAT)



£49
(excl. VAT)

ECG option for Pocket Plus

Enhance the capabilities of your DefiSign Pocket Plus with the ECG-Display Option, allowing you to monitor the patient's heart rhythm. The inclusion of an ECG display is particularly advantageous in professional settings, providing experienced users with valuable insights into the heart rhythm during the resuscitation process. Furthermore, the ECG data can be recorded and monitored after the administration for further analysis and evaluation.

Bags



Every second counts! Don't waste time in an emergency. With this handy **DefiSign defibrillator carry case**, you can transport the DefiSign defib and accessories. The sturdy shoulder strap ensures that the bag hangs securely over your shoulder.

£39 excl. VAT

This **universal defibrillator backpack** is perfectly suited to easily transport your AED and any other medical equipment or accessories, wherever you go. The large front pocket ensures that all AEDs in our range fit into this eye-catching backpack. The front pocket is secured with a clip closure.

£59 excl. VAT



Electrodes

Shop our wide range of **defibrillator electrode pads** for all brands of Automated External Defibrillators. Replace electrode pads before the expiry date to ensure that they work when you need them the most. Most defibrillator electrode pads have a shelf life of around 1.5 - 3 years.



DefiSign Life Electrode Pads
£54 excl. VAT



Schiller Fred Easy adult electrode pads
£35 excl. VAT



Schiller Fred Easyport adult electrode pads
£49 excl. VAT



Heartsine Samaritan Adult PAD PAK
£99 excl. VAT



Philips Heartstart FRx smart II electrode pads
£44 excl. VAT



Philips Heartstart HS1 Adult electrode pads
£49 excl. VAT

Batteries


Every **defibrillator battery** has an expiration date. Please replace the battery in time to ensure a fully functioning device. A beeping sound or flashing light on the AED also indicates that the battery needs to be replaced. Most AED batteries can easily be replaced by yourself. After inserting the new battery the unit will require a self-test (check user manual for guidance). The device is now ready for use; the status indicator will also confirm this.



BatteryPack DefiSign Pocket Plus /
FRED Easyport Plus
£149 excl. VAT



DefiSign Life battery
£149 excl. VAT



Schiller Fred Easyport Lithium Battery
£89 excl. VAT



Schiller Fred Easy Lithium Battery
£129 excl. VAT



Philips Heartstart HS1 and FRX battery
£119 excl. VAT